

Core practical 9: Determine K_a for a weak acid

Objective	
<ul style="list-style-type: none"> To determine K_a for a weak acid 	
Safety	Specification links
<ul style="list-style-type: none"> Wear a lab coat and use eye protection. Tie long hair back. Sodium hydroxide solution is an irritant. 	<ul style="list-style-type: none"> Practical techniques 1, 3, 4, 6, 11 CPAC 1a, 2a, 2b, 3a, 4a, 4b
Procedure	Notes on procedure
<ol style="list-style-type: none"> Set up the datalogger to read the pH, or calibrate the pH meter. Pipette 25.0 cm³ of 0.1 mol dm⁻³ ethanoic acid solution into a 250 cm³ conical flask. Fill a burette with sodium hydroxide solution. Add two or three drops of phenolphthalein to the conical flask. Titrate the ethanoic acid with sodium hydroxide solution until the mixture <i>just</i> turns pink. Pipette a further 25.0 cm³ of 0.1 mol dm⁻³ ethanoic acid solution into the 250 cm³ conical flask. Record the pH of this solution. 	<ul style="list-style-type: none"> Explain to students that titre values should be recorded to 2 decimal places with the second figure being 0 or 5 only. Explain to students that the titration should be repeated until concordant results are obtained.
Answers to questions	
<ol style="list-style-type: none"> Depends on students' findings but the pH should be 4.77; so $[H^+] = 1.7 \times 10^{-5}$ mol dm⁻³ 1.7×10^{-5} mol dm⁻³ Sources of uncertainty include inaccuracy of burette readings, and difficulty identifying the exact end-point. Read glassware from the bottom of the meniscus; use a white tile so you can see the colour change clearly. 	
Sample data	
The pH of the solution will be 4.77, though students' results may vary.	

Core practical 9: Determine K_a for a weak acid**Objective**

- To determine K_a for a weak acid

Safety

- Wear a lab coat and use eye protection.
- Tie long hair back.
- Sodium hydroxide solution is an irritant.

All the maths you need

- Use logarithms in relation to quantities that range over several orders of magnitude.
- Change the subject of an equation.
- Substitute numerical values into algebraic equations using appropriate units for physical quantities.

Equipment

- 100 cm³ of 0.1 mol dm⁻³ ethanoic acid solution
- 100 cm³ of 0.1 mol dm⁻³ sodium hydroxide solution
- datalogger and pH probe or pH meter
- stand, clamp and boss for pH probe
- 50 cm³ burette
- burette stand
- 250 cm³ conical flask
- 25 cm³ pipette and filler
- phenolphthalein indicator

Procedure

1. Set up the datalogger to read the pH, or calibrate the pH meter.
2. Pipette 25.0 cm³ of 0.1 mol dm⁻³ ethanoic acid solution into a 250 cm³ conical flask.
3. Fill a burette with sodium hydroxide solution.
4. Add two or three drops of phenolphthalein to the conical flask.
5. Titrate the ethanoic acid solution with sodium hydroxide solution until the mixture *just* turns pink.
6. Pipette a further 25.0 cm³ of 0.1 mol dm⁻³ ethanoic acid solution into the 250 cm³ conical flask.
7. Record the pH of this solution.

Analysis of results

- Record the pH of the solution after the second addition of ethanoic acid.

Learning tips

- The ionisation of an acid is shown by:



Because there is an equilibrium set up, an equilibrium constant, K_a , can be written:

$$K_a = \frac{[\text{H}^+][\text{A}^-]}{[\text{HA}]}$$

- The K_a value is an indication of acid strength. The larger the value of the K_a , the stronger the acid.
- The K_a of a weak acid can be measured by titrating a known volume of the acid against sodium hydroxide using phenolphthalein as an indicator. A further equal volume of acid is then added, and the pH of the resulting solution is measured. Because effectively half of the acid has been titrated:

$$[\text{H}^+] = [\text{HA}] = [\text{A}^-]$$

$[\text{A}^-]$ and $[\text{HA}]$ can be cancelled in the K_a expression

$$\text{And so } K_a = [\text{H}^+]$$

The pH value of the combined solutions can be converted to $[\text{H}^+]$ to give a K_a value.

Questions

- Use the pH of your solution to calculate $[\text{H}^+]$.
- Calculate a value of K_a for ethanoic acid.
- What are some of the sources of uncertainty in this experiment? What can you do to overcome them?

Exam-style questions

- Folic acid behaves as a weak acid – it can be represented as HA.
Write the equation for the ionisation of a weak acid HA in water. Write the terms *conjugate acid* and *conjugate base* under the relevant formulae for an acid–base pair in your equation. (2)
- Write the expression for the acidity constant K_a for the ionisation of folic acid. (1)
 - If $K_a = 5.0 \times 10^{-3} \text{ mol dm}^{-3}$, calculate $\text{p}K_a$. (1)
 - Calculate the pH of a 0.10 mol dm^{-3} solution of this acid. (2)

Core practical 9: Determine K_a for a weak acid**Objective**

- To determine K_a for a weak acid

Safety

- Wear a lab coat and use eye protection.
- Consult CLEAPSS Hazcards® 38A and 91. Perform a risk assessment using up-to-date information before this practical is carried out.

Equipment per student/group**Notes on equipment**

100 cm³ of 0.1 mol dm⁻³ ethanoic acid solution

Low hazard
If prepared from fresh, remember that pure ethanoic acid is corrosive.

100 cm³ of 0.1 mol dm⁻³ sodium hydroxide solution

Irritant

datalogger and pH probe or pH meter

stand, clamp and boss for pH probe

50 cm³ burette

burette stand

250 cm³ conical flask

25 cm³ pipette and filler

phenolphthalein indicator

Notes